



AMRUT PLUS

Key Features

- ✓ Detailed voice prompts and visual icons for guidance to user
- ✓ Battery capacity : 200 shocks or 4 hours
- ✓ 20 sec delay time between two shocks
- ✓ Evaluates patient ECG to determine if rhythm is shockable
- ✓ 3.5 inch display is for 'Analyzing ', 'Attach Pads', 'CPR', shock status messages, and ECG trace



[Click Here](#) For Technical Specification



*Due to our continuous product development program, the specifications may change without prior notice.

NASAN MEDICAL ELECTRONICS PVT. LTD.

C Wing - Shubham Heights, Off Mumbai - Bangalore Highway Service Road, Warje, Pune - 411058

Phone : +91 8237032389 | E-mail : sales@nasanmedical.com

www.nasanmedical.com