



DEDICA

Dedicated Stress Test System

Key Features

- ✓ 17" colour TFT display with Touch screen operations
- ✓ Offers Rock steady baseline at heart rate more than 210
- ✓ Offers full frequency response from 0.05 to 150 Hz. as stipulated by American Heart Association.
- ✓ Offers accurate computation of ST levels, assuring you of excellent correlation with angiography.
- ✓ Offers clean ECG without muscle tremor or 50 cycle noise.
- ✓ ST integral of all 12 leads

Technical Specifications

Specifications

Length	: 2080 mm.
Width	: 785 mm.
Walking Area	: 1350 X 500 mm.
Weight	: 100 kg.
Belt Speed	: smooth, 0 - 15 kmph.
Incline Range	: 0 - 22 %
Patient Weight	: Upto 250 kg.
Power	: 230V, 50Hz, 15A
In-House Manufacturing	: YES

Features

- ✓ During acquisition when in auto mode, the enlarged median having maximum ST depression is automatically selected and displayed.
- ✓ Enlarged median can be selected manually as well.
- ✓ During acquisition, medians can be refreshed at one second interval from 1 to 8 seconds.
- ✓ Data can be reviewed beat to beat with full disclosure for all 12 leads of the entire test.
- ✓ There is virtually no constraint on the number of cases which can be stored.

Features of Treadmill

- Wide and long belt for stable and comfortable stride.
- 20 mm thick wooden floating deck for longer life.
- Anti-skid side steps.
- Light weight.
- Roller fitted with sealed dust proof ball bearings.
- Strong steel platform design for shock absorption and enhanced cushioning.
- Drive System: 2.0 HP continuous-duty Permanent Magnet DCmotor which ensures trouble-free running even under heavy loads.
- Motor with 4 Kg dynamically balanced flywheel for smooth operation.
- AC motor driven integrated front incline system with built in feedback for correct incline positioning.
- Smooth, noise free operations.



[Click Here For Technical Specification](#)

*Due to our continuous product development program, the specifications may change without prior notice.




12 Lead mode with Superimposed Medians



4 Lead mode with ST Integral



3 Lead mode



NASAN MEDICAL ELECTRONICS PVT. LTD.

C Wing - Shubham Heights, Off Mumbai - Bangalore Highway Service Road, Warje, Pune - 411058

Phone : +91 8237032389 | E-mail : sales@nasanmedical.com

www.nasanmedical.com