

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Linked Median Report

August 19, 2025                    02:02 pm

PATIENT NAME : Demo Patient

Report time                    :            02:12 pm

PROTOCOL        : BRUCE

STAGE            : Pre-Test

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

PHASE TIME        : 00:16

STAGE DURATION   : 00:00

HR                 : 73 bpm

QRS Lead          : V5

THR                :        199

SPEED: 0.0 kmph

GRADE: 0.0

BP                 : 0 / 0 mmHg

METs              : 1.00




I 0.61 / 0.37

aVR -0.68 / -0.67

V1 0.56 / -0.18

V4 1.15 / 0.49




II 0.76 / 0.37

aVL 0.24 / 0.24

V2 3.59 / 2.69

V5 0.56 / 0.24



III 0.12 / -0.06

aVF 0.46 / 0.12

V3 1.86 / 0.85

V6 0.42 / 0.37



V5

BARVE  
MD

10mm/mv, 25mm/sec

100 Hz

NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Linked Median Report

August 19, 2025                    02:04 pm

PATIENT NAME : Demo Patient

Report time                    :            02:12 pm

PROTOCOL        : BRUCE

STAGE            : Exercise Stage 3

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

PHASE TIME        : 00:47

STAGE DURATION : 00:08

HR                 : 74 bpm

QRS Lead          : V5

THR                :        199

SPEED: 5.5 kmph

GRADE: 14.0

BP                 : 0 / 0 mmHg

METs              : 3.40




I 0.20 / -0.31

aVR -0.20 / -0.49

V1 0.44 / 0.06

V4 0.88 / 0.18




II 0.29 / -0.12

aVL 0.22 / 0.12

V2 3.15 / 2.69

V5 0.22 / -0.18




III -0.12 / -0.61

aVF 0.12 / -0.31

V3 1.66 / 0.79

V6 0.05 / -0.24



V5

BARVE  
MD

10mm/mv, 25mm/sec

100 Hz

NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Linked Median Report

August 19, 2025                      02:05 pm

PATIENT NAME : Demo Patient

Report time                      :                      02:12 pm

PROTOCOL        : BRUCE

STAGE            : Peak Exercise

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

PHASE TIME        : 01:47

STAGE DURATION : 00:16

HR                        : 77 bpm

QRS Lead            : V5


THR                      :                      199

SPEED: 9.7 kmph

GRADE: 22.0

BP                        : 0 / 0 mmHg

METs                    : 4.80




I 0.15 / -0.31

aVR -0.07 / -0.24

V1 0.37 / -0.06

V4 0.83 / 0.24




II 0.12 / -0.12

aVL 0.20 / 0.00

V2 3.05 / 2.63

V5 0.20 / -0.06




III -0.22 / -0.61

aVF -0.02 / -0.24

V3 1.61 / 0.79

V6 0.00 / -0.24




V5

BARVE  
MD

10mm/mv, 25mm/sec

100 Hz

NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00



|                                                                                                                  |  |                                                                  |  |                                                                                  |  |
|------------------------------------------------------------------------------------------------------------------|--|------------------------------------------------------------------|--|----------------------------------------------------------------------------------|--|
| PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg                                                                |  | Mixed Median Report                                              |  | August 19, 2025    02:05 pm                                                      |  |
| PATIENT NAME : Demo Patient                                                                                      |  |                                                                  |  | Report time    :    02:12 pm                                                     |  |
| PROTOCOL    : BRUCE<br>STAGE        : Peak Exercise<br>ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J |  | DEMO<br>Heart Rate : 77 (bpm)<br>BP    : 0 / 0 mmHg<br>THR : 199 |  | PHASE TIME    : 01:47<br>STAGE TIME    : 00:16<br>SPEED: 9.7 kmph<br>GRADE: 22.0 |  |

100 Hz

BARVE MD

BARVE  
MD

10mm/mv, 25mm/sec 100 Hz NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00

10mm/mv, 25mm/sec 100 Hz NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00

10mm/mv, 25mm/sec 100 Hz NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Unaveraged ECG Report

August 19, 2025    02:02 pm

PATIENT NAME : Demo Patient


Report time    :    02:12 pm

PROTOCOL    : BRUCE  
STAGE        : Pre-Test

DEMO

PHASE TIME    : 00:16    SPEED: 0.0 kmph  
STAGE DURATION : 00:00    GRADE: 0.0  
HR             : 73 bpm    BP    : 0 / 0 mmHg

THR            :    199



10mm/mv, 25mm/sec

100 Hz

NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Rhythm Report

August 19, 2025                    02:02 pm


PATIENT NAME : Demo Patient

Report time                    :            02:12 pm

PROTOCOL        : BRUCE  
STAGE            : Pre-Test

DEMO

PHASE TIME        : 00:16                    SPEED: 0.0 kmph  
STAGE DURATION : 00:00                    GRADE: 0.0  
HR                 : 73 bpm                    BP        : 0 / 0 mmHg  
QRS Lead          : V5  
THR                :            199




10mm/mv, 25mm/sec

100 Hz

NASAN (C) ST-Win Standard BL 9.0, Unit Ver 0.00







PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Trend Report

August 19, 2025

02:02 pm

PATIENT NAME : Demo Patient

Report time

: 02:12 pm

PROTOCOL : BRUCE

DEMO

I (10mm/mV)




II (10mm/mV)

III (10mm/mV)

ST Level (mm)

ST Level (mm)




ST Level (mm)



ST Slope (mV/sec)

ST Slope (mV/sec)

ST Slope (mV/sec)



PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg  
PATIENT NAME : Demo Patient

Superimpose Report

August 19, 2025    02:02 pm  
Report time        :    02:12 pm

(A) BASE STAGE :  
Wait For Exercise  
00:29 min

77 bpm  
DEMO  
0 / 0

0.00 Kmph

Grade: 0.0 %

PROTOCOL  
BRUCE


(B) COMPARED WITH :

Peak Exercise  
00:19 min

77 bpm  
0 / 0

9.70 Kmph

Grade: 22.0 %



10mm/mv,25mm/sec

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg  
PATIENT NAME : Demo Patient

Superimpose Report


August 19, 2025    02:02 pm  
Report time        :    02:12 pm

(A) BASE STAGE :                      Exercise Stage 1                      72 bpm                      2.70 Kmph                      Grade: 10.0 %                      PROTOCOL

00:18 min                      DEMO                      0 / 0                      BRUCE

(B) COMPARED WITH :                      Peak Exercise                      77 bpm                      9.70 Kmph                      Grade: 22.0 %

00:19 min                      0 / 0



10mm/mv,25mm/sec


02:02 pm

02:12 pm

PROTOCOL  
BRUCE


80 ms Post J

MAX



(0.85)  
(0.24)

MAX



(0.22)  
(-0.12)

MAX



(0.02)  
(-0.24)

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Full Disclosure

August 19, 2025    02:04 pm


PATIENT NAME : Demo Patient

Report time        :    02:12 pm

PROTOCOL        : BRUCE  
STAGE            : Exercise Stage 3

DEMO

PHASE TIME        : 00:47    SPEED: 5.5 kmph  
STAGE DURATION : 00:08    GRADE: 14.0  
THR / HR          :199 / 77 bpm    BP        : 0 / 0 mmHg  
QRS Lead          : V5



10mm/mv, 25mm/sec

100 Hz

DMH

PUNE

PATIENT ID :  
PATIENT NAME : Demo Patient

Summary Report

Report time : 02:12 pm  
August 19, 2025 02:04 pm

DEMO

PROTOCOL : BRUCE

PATIENT HEIGHT : 0 Cm

PATIENT WEIGHT : 0.00 Kg

PATIENT ADD. :

Ref. By : Not Applicable

( Not Applicable )

OBJECT OF TEST :  
RISK FACTOR :  
ACTIVITY :  
MEDICATION :  
BRIEF HISTORY :  
OTHER INVESTIGATION :  
REASON FOR TERMINATION :  
EXERCISE TOLERANCE :  
EXERCISE INDUCED ARRHYTHMIA :  
HAEMO RESPONSE :  
CHRONO RESPONSE :  
FINAL IMPRESSION :

**BARVE  
MD**

DMH

PUNE

PATIENT ID :  
PATIENT NAME : Demo Patient

Summary Report

Report time : 02:12 pm  
August 19, 2025 02:04 pm

DEMO

PROTOCOL : BRUCE

Ref. By : Not Applicable

PATIENT HEIGHT : 0 Cm PATIENT WEIGHT : 0.00 Kg

( Not Applicable )

PATIENT ADD. :

| Stage             | Stage Time | Speed (Kmph) / Grade ( % ) | HR bpm | BP mmHg | R.P.P. X 1000 | METS | ST Level I | Stage Comments |
|-------------------|------------|----------------------------|--------|---------|---------------|------|------------|----------------|
| Pre-Test          | 00:20      | 0.00 / 0.00                | 73     | 0 / 0   | 0             | 1.00 | 0.61       |                |
| Supine            | 00:10      | 0.00 / 0.00                | 77     | 0 / 0   | 0             | 1.00 | 0.54       |                |
| Standing          | 00:10      | 0.00 / 0.00                | 77     | 0 / 0   | 0             | 1.00 | 0.46       |                |
| HyperVentilation  | 00:12      | 0.00 / 0.00                | 74     | 0 / 0   | 0             | 1.00 | 0.46       |                |
| Wait For Exercise | 00:29      | 0.00 / 0.00                | 77     | 0 / 0   | 0             | 1.00 | 0.29       |                |
| Exercise Stage 1  | 00:18      | 2.70 / 10.00               | 73     | 0 / 0   | 0             | 3.40 | 0.29       |                |
| Exercise Stage 2  | 00:18      | 4.00 / 12.00               | 77     | 0 / 0   | 0             | 3.40 | 0.20       |                |
| Exercise Stage 3  | 00:10      | 5.50 / 14.00               | 74     | 0 / 0   | 0             | 3.40 | 0.20       |                |
| Exercise Stage 4  | 00:13      | 6.80 / 16.00               | 73     | 0 / 0   | 0             | 4.80 | 0.20       |                |
| Exercise Stage 5  | 00:14      | 8.10 / 18.00               | 74     | 0 / 0   | 0             | 4.80 | 0.15       |                |
| Exercise Stage 6  | 00:12      | 8.90 / 20.00               | 74     | 0 / 0   | 0             | 4.80 | 0.15       |                |
| Peak Exercise     | 00:19      | 9.70 / 22.00               | 77     | 0 / 0   | 0             | 4.80 | 0.15       |                |
| Recovery          | 00:22      | 0.00 / 0.00                | 74     | 0 / 0   | 0             | 1.00 | 0.15       |                |

TOTAL EXER TIME : 1 : 44 min

MAX HR : 77 bpm (38.69 % of 199 bpm)

MAX BP : 0 / 0 mmHg

MAX WORKLOAD : 4.80

DISTANCE COVERED : 0.19 Km

DOUBLE PRODUCT : 0.00

**BARVE**  
**MD**

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

ST Tables

August 19, 2025

02:04 pm

PATIENT NAME : Demo Patient

Report time        :

02:12 pm

PROTOCOL        : BRUCE

## DEMO

## ST LEVELS

| Phase             | I   | II  | III  | aVR  | aVL | aVF  | V1  | V2  | V3  | V4  | V5  | V6  |
|-------------------|-----|-----|------|------|-----|------|-----|-----|-----|-----|-----|-----|
| Pre-Test          | 0.6 | 0.8 | 0.1  | -0.7 | 0.2 | 0.5  | 0.6 | 3.6 | 1.9 | 1.1 | 0.6 | 0.4 |
| Supine            | 0.5 | 0.7 | 0.1  | -0.6 | 0.2 | 0.4  | 0.6 | 3.5 | 1.7 | 1.0 | 0.5 | 0.4 |
| Standing          | 0.5 | 0.6 | 0.0  | -0.5 | 0.2 | 0.3  | 0.5 | 3.3 | 1.6 | 0.9 | 0.4 | 0.3 |
| HyperVentilation  | 0.5 | 0.6 | 0.0  | -0.5 | 0.2 | 0.3  | 0.5 | 3.3 | 1.7 | 1.0 | 0.4 | 0.3 |
| Wait For Exercise | 0.3 | 0.4 | -0.0 | -0.3 | 0.2 | 0.2  | 0.5 | 3.2 | 1.6 | 0.9 | 0.2 | 0.2 |
| Exercise Stage 1  | 0.3 | 0.4 | -0.0 | -0.3 | 0.2 | 0.2  | 0.5 | 3.2 | 1.7 | 0.9 | 0.2 | 0.2 |
| Exercise Stage 2  | 0.2 | 0.3 | -0.1 | -0.2 | 0.2 | 0.1  | 0.4 | 3.1 | 1.6 | 0.9 | 0.2 | 0.0 |
| Exercise Stage 3  | 0.2 | 0.3 | -0.1 | -0.2 | 0.2 | 0.1  | 0.4 | 3.2 | 1.7 | 0.9 | 0.2 | 0.0 |
| Exercise Stage 4  | 0.2 | 0.3 | -0.1 | -0.2 | 0.2 | 0.1  | 0.4 | 3.2 | 1.6 | 0.8 | 0.2 | 0.0 |
| Exercise Stage 5  | 0.1 | 0.2 | -0.2 | -0.1 | 0.2 | 0.0  | 0.4 | 3.1 | 1.6 | 0.8 | 0.2 | 0.0 |
| Exercise Stage 6  | 0.1 | 0.2 | -0.2 | -0.1 | 0.2 | 0.0  | 0.4 | 3.1 | 1.7 | 0.9 | 0.2 | 0.0 |
| Peak Exercise     | 0.1 | 0.1 | -0.2 | -0.1 | 0.2 | -0.0 | 0.4 | 3.1 | 1.6 | 0.8 | 0.2 | 0.0 |
| Recovery          | 0.1 | 0.1 | -0.2 | -0.1 | 0.2 | 0.0  | 0.4 | 3.1 | 1.7 | 0.9 | 0.2 | 0.0 |

## ST SLOPES

| Phase             | I    | II   | III  | aVR  | aVL | aVF  | V1   | V2  | V3  | V4  | V5   | V6   |
|-------------------|------|------|------|------|-----|------|------|-----|-----|-----|------|------|
| Pre-Test          | 0.4  | 0.4  | -0.1 | -0.7 | 0.2 | 0.1  | -0.2 | 2.7 | 0.9 | 0.5 | 0.2  | 0.4  |
| Supine            | 0.4  | 0.4  | -0.1 | -0.6 | 0.2 | 0.1  | -0.2 | 2.7 | 0.9 | 0.5 | 0.2  | 0.4  |
| Standing          | 0.1  | 0.1  | -0.3 | -0.7 | 0.2 | -0.1 | -0.1 | 2.7 | 0.7 | 0.4 | 0.1  | 0.1  |
| HyperVentilation  | 0.1  | 0.1  | -0.3 | -0.7 | 0.2 | -0.1 | -0.1 | 2.7 | 0.7 | 0.4 | 0.1  | 0.1  |
| Wait For Exercise | -0.2 | 0.0  | -0.5 | -0.5 | 0.2 | -0.2 | 0.0  | 2.6 | 0.7 | 0.1 | -0.1 | -0.1 |
| Exercise Stage 1  | -0.2 | 0.1  | -0.4 | -0.4 | 0.1 | -0.1 | -0.1 | 2.7 | 0.8 | 0.2 | 0.0  | 0.0  |
| Exercise Stage 2  | -0.3 | -0.1 | -0.6 | -0.5 | 0.1 | -0.3 | 0.1  | 2.6 | 0.7 | 0.1 | -0.2 | -0.2 |
| Exercise Stage 3  | -0.3 | -0.1 | -0.6 | -0.5 | 0.1 | -0.3 | 0.1  | 2.7 | 0.8 | 0.2 | -0.2 | -0.2 |
| Exercise Stage 4  | -0.3 | 0.0  | -0.5 | -0.3 | 0.1 | -0.2 | 0.0  | 2.7 | 0.9 | 0.2 | -0.1 | -0.2 |
| Exercise Stage 5  | -0.5 | -0.2 | -0.7 | -0.4 | 0.1 | -0.4 | -0.2 | 2.4 | 0.5 | 0.0 | -0.2 | -0.4 |
| Exercise Stage 6  | -0.4 | -0.2 | -0.7 | -0.4 | 0.1 | -0.4 | -0.1 | 2.7 | 0.8 | 0.2 | -0.2 | -0.3 |
| Peak Exercise     | -0.3 | -0.1 | -0.6 | -0.2 | 0.0 | -0.2 | -0.1 | 2.6 | 0.8 | 0.2 | -0.1 | -0.2 |
| Recovery          | -0.3 | -0.1 | -0.6 | -0.2 | 0.0 | -0.2 | -0.1 | 2.7 | 0.8 | 0.2 | -0.1 | -0.2 |

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Extended Mixed Median Report

August 19, 2025    02:04 pm

PATIENT NAME : Demo Patient

Report time        :    02:12 pm

PROTOCOL        : BRUCE  
STAGE            : Pre-Test  
ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

Heart Rate : 73 (bpm)

PHASE TIME        : 00:16

SPEED: 0.0 kmph

THR                : 199

STAGE TIME        : 00:16

GRADE: 0.0

BP                 : 0 / 0 mmHg



BARVE  
MD

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

Extended Mixed Median Report

August 19, 2025    02:04 pm

PATIENT NAME : Demo Patient

Report time    :    02:12 pm

PROTOCOL    : BRUCE  
STAGE        : Exercise Stage 3  
ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

Heart Rate : 74 (bpm)

PHASE TIME    : 00:47

SPEED: 5.5 kmph

THR    : 199

STAGE TIME    : 00:08

GRADE: 14.0

BP    : 0 / 0 mmHg



BARVE  
MD

PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg


Selected Median Report

August 19, 2025    02:04 pm

PATIENT NAME : Demo Patient

Report time    :    02:12 pm

Protocol : BRUCE    Total Exercise Time : 1 : 44 min    Max HR : 77 bpm (38.69 % of 199 bpm)    Max. BP : 0 / 0 mm/Hg    Max. Workload : 4.80  
DEMO



PATIENT ID :    HEIGHT : 0 Cm    WEIGHT : 0.00 Kg

MAX HR Report

August 19, 2025                    02:04 pm


PATIENT NAME : Demo Patient

Report time                    :            02:12 pm

PROTOCOL        : BRUCE  
STAGE            : Pre-Test

DEMO

PHASE TIME        : 00:08                    SPEED: 0.0 kmph  
STAGE DURATION : 00:00                    GRADE: 0.0  
THR / HR           :199 / 77 bpm            BP        : 0 / 0 mmHg  
QRS Lead           : V5



10mm/mv,25mm/sec

100 Hz