

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Linked Median Report

October 13, 2024 03:03 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : Pre-Test

DEMO

PHASE TIME : 00:08

SPEED: 0.0 kmph

STAGE DURATION : 00:00

GRADE: 0.0

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J


HR : 76 bpm

BP : 0 / 0 mmHg

QRS Lead : V5

METs : 1.00

THR : 164




I 0.63 / 0.37

aVR -0.71 / -0.85

V1 0.59 / -0.12

V4 1.29 / 0.43




II 0.78 / 0.24

aVL 0.24 / 0.31

V2 3.61 / 2.69

V5 0.61 / 0.18




III 0.12 / -0.12

aVF 0.46 / 0.06

V3 2.00 / 0.92

V6 0.44 / 0.31



V5

Dr.AAAA
DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Linked Median Report

October 13, 2024 03:03 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : Supine

DEMO

PHASE TIME : 00:16

SPEED: 0.0 kmph

STAGE DURATION : 00:00

GRADE: 0.0

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J


HR : 74 bpm

BP : 0 / 0 mmHg

QRS Lead : V5

METs : 1.00

THR : 164




I 0.63 / 0.37

aVR -0.71 / -0.85

V1 0.59 / -0.12

V4 1.22 / 0.43




II 0.78 / 0.24

aVL 0.24 / 0.31

V2 3.61 / 2.69

V5 0.59 / 0.18




III 0.12 / -0.12

aVF 0.46 / 0.06

V3 1.93 / 0.92

V6 0.44 / 0.31



V5

Dr.AAAA
DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

October 13, 2024

03:03 pm

PATIENT NAME : Demo Patient

Linked Median Report

Report time

: 03:11 pm

PROTOCOL : BRUCE

STAGE : Standing

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

PHASE TIME : 00:24

SPEED: 0.0 kmph

STAGE DURATION : 00:00

GRADE: 0.0

HR : 73 bpm

BP : 0 / 0 mmHg

QRS Lead : V5

METs : 1.00

THR : 164




I 0.61 / 0.37

aVR -0.68 / -0.67

V1 0.56 / -0.18

V4 1.10 / 0.49




II 0.76 / 0.37

aVL 0.24 / 0.24

V2 3.59 / 2.69

V5 0.56 / 0.24




III 0.12 / -0.06

aVF 0.46 / 0.12

V3 1.81 / 0.85

V6 0.42 / 0.37



V5

Dr.AAAA

DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Linked Median Report

October 13, 2024 03:04 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : HyperVentilation

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

PHASE TIME : 00:40

STAGE DURATION : 00:00

HR : 72 bpm

QRS Lead : V5


THR : 164

SPEED: 0.0 kmph

GRADE: 0.0

BP : 0 / 0 mmHg

METs : 1.00




I 0.56 / 0.12

aVR -0.61 / -0.73

V1 0.56 / -0.06

V4 0.95 / 0.43




II 0.68 / 0.12

aVL 0.27 / 0.24

V2 3.54 / 2.69

V5 0.44 / 0.12




III 0.07 / -0.31

aVF 0.42 / -0.12

V3 1.73 / 0.73

V6 0.34 / 0.12



V5

Dr.AAAA

DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Linked Median Report

October 13, 2024

03:04 pm

PATIENT NAME : Demo Patient

Report time

: 03:11 pm

PROTOCOL : BRUCE

STAGE : Wait For Exercise

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

PHASE TIME : 00:53

STAGE DURATION : 00:00

HR : 72 bpm

QRS Lead : V5

THR : 164

SPEED: 0.0 kmph

GRADE: 0.0

BP : 0 / 0 mmHg

METs : 1.00



I 0.56 / 0.12

aVR -0.61 / -0.73

V1 0.56 / -0.06

V4 0.95 / 0.43



II 0.68 / 0.12

aVL 0.27 / 0.24

V2 3.54 / 2.69

V5 0.44 / 0.12



III 0.07 / -0.31

aVF 0.42 / -0.12

V3 1.76 / 0.79

V6 0.34 / 0.12



V5

Dr.AAAA
DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Linked Median Report

October 13, 2024 03:07 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : Exercise Stage 1

DEMO

PHASE TIME : 03:00

SPEED: 2.7 kmph

STAGE DURATION : 03:04

GRADE: 10.0

HR : 74 bpm


BP : 0 / 0 mmHg

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

QRS Lead : V5

METs : 5.10

THR : 164



I 0.24 / -0.31

aVR -0.22 / -0.49

V1 0.46 / 0.06

V4 0.85 / 0.12




II 0.32 / -0.12

aVL 0.24 / 0.06

V2 3.20 / 2.56

V5 0.20 / -0.06




III -0.10 / -0.61

aVF 0.12 / -0.37

V3 1.66 / 0.73

V6 0.05 / -0.24



V5

Dr.AAAA
DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Linked Median Report

October 13, 2024 03:07 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : Peak Exercise

DEMO

PHASE TIME : 03:23

SPEED: 4.0 kmph

STAGE DURATION : 00:16

GRADE: 12.0

HR : 75 bpm




BP : 0 / 0 mmHg

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

QRS Lead : V5

METs : 5.70

THR : 164



Dr.AAAA
DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Linked Median Report

October 13, 2024 03:08 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : Recovery 1

DEMO

PHASE TIME : 04:23

SPEED: 1.2 kmph

STAGE DURATION : 01:04

GRADE: 0.0

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J




HR : 74 bpm

BP : 123 / 78 mmHg

QRS Lead : V5

METs : 6.60

THR : 164



Comments :

Dr.AAAA
DM

10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

October 13, 2024 03:08 pm

PATIENT NAME : Demo Patient

Mixed Median Report

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : Recovery 1

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

Heart Rate : 74 (bpm)

BP : 123 / 78 mmHg

THR : 164

PHASE TIME : 04:23

STAGE TIME : 01:04

SPEED: 1.2 kmph

GRADE: 0.0



Dr.AAAA
DM

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Unaveraged ECG Report

October 13, 2024 03:04 pm

PATIENT NAME : Demo Patient


Report time : 03:11 pm

PROTOCOL : BRUCE
STAGE : Wait For Exercise

DEMO

PHASE TIME : 00:53 SPEED: 0.0 kmph
STAGE DURATION : 00:00 GRADE: 0.0
HR : 72 bpm BP : 0 / 0 mmHg

THR : 164



10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

PATIENT NAME : Demo Patient

Rhythm Report


October 13, 2024 03:04 pm

Report time : 03:11 pm

PROTOCOL : BRUCE
STAGE : Wait For Exercise

DEMO

PHASE TIME : 00:53 SPEED: 0.0 kmph
STAGE DURATION : 00:00 GRADE: 0.0
HR : 72 bpm BP : 0 / 0 mmHg
QRS Lead : V5
THR : 164



10mm/mv, 25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00















PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

PATIENT NAME : Demo Patient

Stage Report

October 13, 2024 03:04 pm

Report time : 03:11 pm

PROTOCOL	: BRUCE							
	Pre-Test	Supine	Standing	DEMO	HyperVentilation	Wait For Exercise	Exercise Stage 1	Peak Exercise
Clock1/Clock2	: 00:11/00:00	00:21/00:00	00:32/00:00	00:43/00:00	00:53/00:00	03:00/03:00	03:23/00:22	
Speed/Grade	: 0.00 Kmph/0.00%	0.00 Kmph/0.00%	0.00 Kmph/0.00%	0.00 Kmph/0.00%	0.00 Kmph/0.00%	2.70 Kmph/10.00%	4.00 Kmph/12.00%	
HR / BP	: 76 bpm 0/0	74 bpm 0/0	73 bpm 0/0	72 bpm 0/0	72 bpm 0/0	74 bpm 0/0	75 bpm 0/0	
ST Levels (mm)								
	V5 (0.61)	V5 (0.59)	V5 (0.56)	V5 (0.44)	V5 (0.44)	V5 (0.20)	V5 (0.20)	
								
	V5 (0.61)	V5 (0.59)	V5 (0.56)	V5 (0.44)	V5 (0.44)	V5 (0.20)	V5 (0.20)	
								
	V5 (0.61)	V5 (0.59)	V5 (0.56)	V5 (0.44)	V5 (0.44)	V5 (0.20)	V5 (0.20)	
								
	V5 (0.61)	V5 (0.59)	V5 (0.56)	V5 (0.44)	V5 (0.44)	V5 (0.20)	V5 (0.20)	
								
	V5 (0.61)	V5 (0.59)	V5 (0.56)	V5 (0.44)	V5 (0.44)	V5 (0.20)	V5 (0.20)	
								
	V5 (0.61)	V5 (0.59)	V5 (0.56)	V5 (0.44)	V5 (0.44)	V5 (0.20)	V5 (0.20)	
						10mm/mv, 25mm/sec		

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Stage Report

October 13, 2024

03:04 pm

PATIENT NAME : Demo Patient

Report time :

03:11 pm

PROTOCOL : BRUCE

Recovery 1

DEMO

Clock1/Clock2 : 04:23/01:00

Speed/Grade : 1.20 Km/h/0.00%

HR / BP : 74 bpm 123/78

ST Levels (mm)

V5 (0.24)

V5 (0.24)

V5 (0.24)

V5 (0.24)

V5 (0.24)

V5 (0.24)

10mm/mv,25mm/sec

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

PATIENT NAME : Demo Patient

Trend Report

October 13, 2024

03:04 pm

Report time :

03:11 pm

PROTOCOL : BRUCE

DEMO

I (10mm/mV)






II (10mm/mV)

III (10mm/mV)

J Amp (mm)

J Amp (mm)

J Amp (mm)



PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

October 13, 2024

03:04 pm

PATIENT NAME : Demo Patient

Trend Report

Report time :

03:11 pm

PROTOCOL : BRUCE

DEMO

I (10mm/mV)




II (10mm/mV)

III (10mm/mV)

ST Level (mm)

ST Level (mm)




ST Level (mm)



ST Slope (mV/sec)

ST Slope (mV/sec)

ST Slope (mV/sec)



PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

October 13, 2024

03:04 pm

PATIENT NAME : Demo Patient

Trend Report

Report time :

03:11 pm

PROTOCOL : BRUCE

DEMO

I (10mm/mV)




II (10mm/mV)

III (10mm/mV)


R Amp (mm)

R Amp (mm)


R Amp (mm)




PVC (count / min)



METS



RPP



PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Superimpose Report

October 13, 2024

03:04 pm

PATIENT NAME : Demo Patient

Report time :

03:11 pm

(A) BASE STAGE :

Wait For Exercise
00:10 min

72 bpm
DEMO
0 / 0

0.00 Kmph

Grade: 0.0 %

PROTOCOL
BRUCE









(B) COMPARED WITH :

Exercise Stage 1
03:00 min

72 bpm
0 / 0

2.70 Kmph

Grade: 10.0 %



10mm/mv,25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg
PATIENT NAME : Demo Patient

Brief Summary Report


October 13, 2024 03:04 pm
Report time : 03:11 pm

BASE STAGE : Supine 74 bpm 0.00 METS 80 ms Post J PROTOCOL
00:10 DEMO 0/0 BRUCE

MAX STAGE : Recovery 1 72 bpm 1.00 METS 80 ms Post J

01:00

123 / 78



10mm/mv,25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

October 13, 2024 03:06 pm

PATIENT NAME : Demo Patient


Full Disclosure

Report time : 03:11 pm

PROTOCOL : BRUCE
STAGE : Exercise Stage 1

DEMO

PHASE TIME : 03:00 SPEED: 2.7 kmph
STAGE DURATION : 02:24 GRADE: 10.0
THR / HR :164 / 75 bpm BP : 0 / 0 mmHg
QRS Lead : V5



10mm/mv,25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID :
PATIENT NAME : Demo Patient

Summary Report

Report time : 03:11 pm
October 13, 2024 03:06 pm

DEMO

PROTOCOL : BRUCE

PATIENT HEIGHT : 0 Cm

PATIENT WEIGHT : 0.00 Kg

PATIENT ADD. :

Ref. By : Not Applicable

(Not Applicable)

OBJECT OF TEST :

RISK FACTOR :

ACTIVITY :

MEDICATION :

BRIEF HISTORY :

OTHER INVESTIGATION :

REASON FOR TERMINATION :

EXERCISE TOLERANCE :

EXERCISE INDUCED ARRHYTHMIA :

HAEMO RESPONSE :

CHRONO RESPONSE :

FINAL IMPRESSION :

Dr.AAAA
DM

PATIENT ID :
PATIENT NAME : Demo Patient

Summary Report

Report time : 03:11 pm
October 13, 2024 03:06 pm

DEMO

PROTOCOL : BRUCE

Ref. By : Not Applicable

PATIENT HEIGHT : 0 Cm

PATIENT WEIGHT : 0.00 Kg

(Not Applicable)

PATIENT ADD. :

Stage	Stage Time	Speed (Kmph) / Grade (%)	HR bpm	BP mmHg	R.P.P. X 1000	METS	ST Level I	Stage Comments
Pre-Test	00:10	0.00 / 0.00	76	0 / 0	0	1.00	0.63	
Supine	00:10	0.00 / 0.00	74	0 / 0	0	1.00	0.63	
Standing	00:10	0.00 / 0.00	73	0 / 0	0	1.00	0.61	
HyperVentilation	00:10	0.00 / 0.00	72	0 / 0	0	1.00	0.56	
Wait For Exercise	00:10	0.00 / 0.00	72	0 / 0	0	1.00	0.56	
Exercise Stage 1	03:00	2.70 / 10.00	74	0 / 0	0	5.70	0.24	
Peak Exercise	00:22	4.00 / 12.00	75	0 / 0	0	5.70	0.22	
Recovery 1	01:00	1.20 / 0.00	74	123 / 78	9	6.60	0.17	

TOTAL EXER TIME : 4 : 22 min

MAX HR : 76 bpm (46.34 % of 164 bpm)

MAX BP : 123 / 78 mmHg

MAX WORKLOAD : 6.60

DISTANCE COVERED : 0.18 Km

DOUBLE PRODUCT : 9348.00

Dr.AAAA
DM

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

ST Tables

October 13, 2024

03:06 pm

PATIENT NAME : Demo Patient

Report time :

03:11 pm

PROTOCOL : BRUCE

DEMO

ST LEVELS

Phase	I	II	III	aVR	aVL	aVF	V1	V2	V3	V4	V5	V6
Pre-Test	0.6	0.8	0.1	-0.7	0.2	0.5	0.6	3.6	2.0	1.3	0.6	0.4
Supine	0.6	0.8	0.1	-0.7	0.2	0.5	0.6	3.6	1.9	1.2	0.6	0.4
Standing	0.6	0.8	0.1	-0.7	0.2	0.5	0.6	3.6	1.8	1.1	0.6	0.4
HyperVentilation	0.6	0.7	0.1	-0.6	0.3	0.4	0.6	3.5	1.7	1.0	0.4	0.3
Wait For Exercise	0.6	0.7	0.1	-0.6	0.3	0.4	0.6	3.5	1.8	1.0	0.4	0.3
Exercise Stage 1	0.2	0.3	-0.1	-0.2	0.2	0.1	0.5	3.2	1.7	0.9	0.2	0.0
Peak Exercise	0.2	0.3	-0.1	-0.2	0.2	0.1	0.4	3.2	1.6	0.9	0.2	0.0
Recovery 1	0.2	0.2	-0.2	-0.1	0.1	0.0	0.4	3.1	1.7	0.9	0.2	0.1

ST SLOPES

Phase	I	II	III	aVR	aVL	aVF	V1	V2	V3	V4	V5	V6
Pre-Test	0.4	0.2	-0.1	-0.9	0.3	0.1	-0.1	2.7	0.9	0.4	0.2	0.3
Supine	0.4	0.2	-0.1	-0.9	0.3	0.1	-0.1	2.7	0.9	0.4	0.2	0.3
Standing	0.4	0.4	-0.1	-0.7	0.2	0.1	-0.2	2.7	0.9	0.5	0.2	0.4
HyperVentilation	0.1	0.1	-0.3	-0.7	0.2	-0.1	-0.1	2.7	0.7	0.4	0.1	0.1
Wait For Exercise	0.1	0.1	-0.3	-0.7	0.2	-0.1	-0.1	2.7	0.8	0.4	0.1	0.1
Exercise Stage 1	-0.3	-0.1	-0.6	-0.5	0.1	-0.4	0.1	2.6	0.7	0.1	-0.1	-0.2
Peak Exercise	-0.2	-0.1	-0.6	-0.5	0.1	-0.3	0.1	2.7	0.8	0.2	-0.1	-0.2
Recovery 1	-0.4	-0.2	-0.7	-0.4	0.0	-0.4	-0.1	2.7	0.7	0.1	-0.2	-0.3

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Extra Comments

October 13, 2024

03:06 pm

PATIENT NAME : Demo Patient

Report time :

03:11 pm

PROTOCOL : BRUCE

DEMO

Extra Comments :

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Advanced Mixed Median Report

October 13, 2024 03:06 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE
STAGE : Exercise Stage 1

ST level in mm @ J+ 80 ms DEMO

Heart Rate : 74 (bpm)
BP : 0 / 0 mmHg
THR : 164

PHASE TIME : 03:00
STAGE TIME : 03:04

SPEED: 2.7 kmph
GRADE: 10.0



Dr.AAAA
DM

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

October 13, 2024

03:06 pm

PATIENT NAME : Demo Patient


Advanced Mixed Median Report

Report time

: 03:11 pm

PROTOCOL : BRUCE
STAGE : Peak Exercise

ST level in mm @ J+ 80 ms DEMO

Heart Rate : 75 (bpm)
BP : 0 / 0 mmHg
THR : 164PHASE TIME : 03:23
STAGE TIME : 00:16SPEED: 4.0 kmph
GRADE: 12.0

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

Extended Mixed Median Report

October 13, 2024 03:06 pm

PATIENT NAME : Demo Patient

Report time : 03:11 pm

PROTOCOL : BRUCE

STAGE : Exercise Stage 1

ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

Heart Rate : 74 (bpm)

THR : 164


BP : 0 / 0 mmHg

PHASE TIME : 03:00

STAGE TIME : 03:04

SPEED: 2.7 kmph

GRADE: 10.0



Dr.AAAA
DM

PROTOCOL : BRUCE
STAGE : Peak Exercise
ST Levels (mm) / ST Slope (mV/s) measured at 80 ms Post J

DEMO

Heart Rate : 75 (bpm)
THR : 164
BP : 0 / 0 mmHg

PHASE TIME : 03:23
STAGE TIME : 00:16

SPEED: 4.0 kmph
GRADE: 12.0



PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

October 13, 2024

03:06 pm

PATIENT NAME : Demo Patient

Selected Median Report

Report time :

03:12 pm

Protocol : BRUCE

Total Exercise Time : 4 : 22 min

Max HR : 76 bpm (46.34 % of 164 bpm)
DEMO

Max. BP : 123 / 78 mm/Hg

Max. Workload : 6.60

SUPINE

MAX-ST

PEAK

END OF
RECOVERYSUPINE
EXERCISEMAX-ST
EXERCISEPEAK
EXERCISEEND OF
RECOVERY

00:10

01:00

01:00

01:00

00:10

01:00

01:00

01:00

74 bpm

72 bpm

74 bpm

74 bpm

74 bpm

72 bpm

74 bpm

74 bpm

0 / 0

123 / 78

123 / 78



123 / 78

0 / 0

123 / 78

123 / 78

123 / 78



Dr.Name: Dr.AAAA(DM)

10mm/mv,25mm/sec

100 Hz

NASAN (C) DEDICA BL 9.0, Unit Ver 0.00

PATIENT ID : HEIGHT : 0 Cm WEIGHT : 0.00 Kg

MAX HR Report

October 13, 2024 03:06 pm


PATIENT NAME : Demo Patient

Report time : 03:12 pm

PROTOCOL : BRUCE
STAGE : Pre-Test

DEMO

PHASE TIME : 00:16 SPEED: 0.0 kmph
STAGE DURATION : 00:00 GRADE: 0.0
THR / HR :164 / 76 bpm BP : 0 / 0 mmHg
QRS Lead : V5



10mm/mv, 25mm/sec